


**OnSite™**  
**H. Pylori Ag Rapid Test**

*A lateral flow immunoassay for the qualitative detection of Helicobacter pylori  
(H. pylori) Antigen in human fecal specimens*



**CTK Biotech, Inc.**  
13855 Stowe Dr,  
Poway, CA 92064, USA  
Tel: 858-457-8698  
Fax: 858-535-1739  
E-mail: [info@ctkbiotech.com](mailto:info@ctkbiotech.com)



MDSS GmbH,  
Schiffgraben 41.  
30175 Hannover, Germany

### Table of Contents

| Section:   | <b>Page(s)</b> |
|--|----------------|
| A Product Description  | 3              |
| 1. Brief Product Description                                 | 3              |
| 2. Intended Use and Operation                                | 3              |
| 3. Device Classification                                     | 3              |
| B. Technical Requirements & Essential Requirements Checklist | 4              |
| 1. General Requirements                                      | 4              |
| 2. Design and Manufacturing Requirements                     | 5              |
| 3. Risk Analysis   | 6              |
| 4. Labels & Instructions for Use                             | 6              |
| 4.1 Instructions for Use                                     | 6              |
| 4.2 Direct Labeling  | 6              |
| 5. Materials Specifications                                  | 7              |
| 5.1 Material Specifications                                  | 7              |
| 5.2 Product Specifications                                   | 7              |
| 5.3 Component Specifications                                 | 7              |
| 6. Clinical Studies & Performance Evaluations                | 9              |
| 7. Stability Studies   | 9              |
| 8. Manufacturing   | 9              |
| 8.1 Manufacturing Process                                    | 9              |
| 8.2 Manufacturing Flow Chart                                 | 11             |
| 8.3 Production Specifications                                | 12             |
| 8.4 Documentation of Quality System                          | 12             |
| 9. Batch Release Criteria                                    | 13             |
| 9.1 Incoming Material Inspection                             | 13             |
| 9.2 In Process Inspection                                    | 13             |
| 9.3 Final Acceptance Inspection                              | 13             |
| 9.4 Procedure for Final Acceptance Inspection                | 14             |
| 10. Conclusion   | 16             |

**Section A) PRODUCT DESCRIPTION****1. Product Description**

The *OnSite H. Pylori Ag Rapid Test* is a lateral flow chromatographic immunoassay for the qualitative detection of *H. pylori* antigen in human fecal specimen. It is intended to be used by professionals as a screening test and provides a preliminary test result to aid in the diagnosis of infection with *H. pylori*. Any interpretation or use of this preliminary test result must also rely on other clinical findings as well as on the professional judgment of health care providers. Alternative test method(s) should be considered to confirm the test result obtained by this device.

**2. Intended Use**

The *OnSite H. Pylori Ag Rapid Test* is intended to be used as a screening test and provides a preliminary test result to aid in the diagnosis of infection with *H. pylori*.

**3. Device Classification**

- USA: Class II device
- Europe: Annex III Other
- India: Non-critical device

### Section B) TECHNICAL REQUIREMENTS

#### 1. General Requirements Checklist

| General requirements |   | Apply | Applied Standards   | Demonstrated By:                  | Location:                       |
|----------------------|---|-------|---------------------|-----------------------------------|---------------------------------|
| 1.                   | Safe use for Patient, User, Environments  | Yes   | EN ISO 14971 : 2012 | Risk Management File              | Attachment # 1                  |
| 2.                   | Solutions To Ensure Safety, Including Elimination/ reduction of risk, Taking Appropriate Action when appropriate and informing users of residual risk                     | Yes   | EN ISO 14971 : 2012 | Risk Management File              | Attachment # 1                  |
| 3                    | They are suitable for the purposes referred to in Article 1(2)(b), & meet manufacturer stated Performance Expectations  | Yes   | EN 13612 : 2002     | IFU & Product Performance         | Attachment # 2 & Attachment # 4 |
| 4.                   | Product Safety and Performance must not be affected during product lifetime when exposed to normal stresses and conditions  | Yes   | EN ISO 23640 : 2013 | Stability Study                   | Attachment #5                   |
| 5.                   | Devices designed and manufactured so that performance is not adversely affected under storage and transport conditions.   | Yes   | EN ISO 23640 : 2013 | Pouch Package Study               | Attachment # 9                  |
| 6.                   | In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use (ISO 18113-2:2009) | Yes   | EN ISO 18113-2:2009 | IFU                               | Attachment # 2                  |
| 7.                   | In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements                          | Yes   | EN ISO 18113:2009   | IFU                               | Attachment # 2                  |
| 8.                   | Sampling procedures used for acceptance testing of in vitro diagnostic medical devices - Statistical aspects  | Yes   | EN ISO 13975: 2003  | Sampling and acceptance procedure | SOP-82-04 and WI 82-04-1        |

### 2. Design and Manufacturing Requirements

| Design and Manufacturing Requirements |  | Apply     | Standards   | Demonstrated By:                   | Location:                      |
|---------------------------------------|--|-----------|---|------------------------------------|--------------------------------|
| <b>1</b>                              | <b>CHEMICAL AND PHYSICAL PROPERTIES</b>                              |           |   |                                    |                                |
| 1.1                                   | Device characteristics and performance in relation with intended use | Yes       | EN 13612 : 2002   | Product Performance                | Attachment #4                  |
| 1.2                                   | Risks posed by device leakage, contaminants and residues             | Yes       | EN ISO 14971 : 2012<br>EN 13641: 2002                       | MSDS                               | Attachment # 8                 |
| <b>2</b>                              | <b>INFECTION AND MICROBIAL CONTAMINATION</b>                         |           |   |                                    |                                |
| <b>2.1</b>                            | Reduction of the risk of infection and/or contamination              | Yes       | EN ISO 14971 : 2012<br>EN 13641:2002                        | Risk Management, MSDS, IFU         | Attachment #1, #2 & # 8        |
| <b>2.2</b>                            | Reduce risk of biological substances                                 | Yes       | EN ISO 14971 : 2012<br>EN 13641:2002                        | Risk Management, MSDS, IFU         | Attachment #1, #2 & # 8        |
| <b>2.3-2.6</b>                        | N/A  | No        | N/A   | N/A                                | N/A                            |
| <b>2.7</b>                            | Packaging  | Yes       | EN ISO 14971 : 2012<br>EN 13640:2002                        | Pouch Package Study                | Attachment # 9                 |
| <b>3</b>                              | <b>MANUFACTURING AND ENVIRONMENTAL PROPERTIES</b>                    |           |   |                                    |                                |
| <b>3.1</b>                            | Tests when connected with other devices and/or accessories           | No        | N/A   | N/A                                | N/A                            |
| <b>3.2</b>                            | Contact with device materials  | No        | N/A   | N/A                                | N/A                            |
| <b>3.3</b>                            | Remove risk due to outside influence. (ex. Humidity, temperature)    | Yes       | EN ISO 14971 : 2012<br>EN 13640:2002                        | Risk Management<br>Stability study | Attachment # 1 & Attachment #5 |
| <b>3.4</b>                            | Flammability   | No        | N/A   | N/A                                | N/A                            |
| <b>3.5</b>                            | Safe Waste Disposal  | Yes       | EN ISO 14971 : 2012<br>EN 13640:2002                        | IFU, MSDS                          | Attachment # 2 & Attachment #8 |
| <b>3.6</b>                            | Ergonomics   | No        | N/A   | N/A                                | N/A                            |
| <b>4-7</b>                            |  | <b>No</b> | <b>N/A</b>  | <b>N/A</b>                         | <b>N/A</b>                     |
| <b>8</b>                              | <b>MANUFACTURER INFORMATION</b>                                      |           |   |                                    |                                |
| <b>8.1</b>                            | Information for use  | Yes       | EN1041 : 2008<br>EN ISO 18113 : 2011<br>EN ISO 15223-1:2012 | Labels, IFU                        | Section 4.2 & Attachment #2    |
| <b>8.2</b>                            | Standard Symbol  | Yes       | EN ISO 15223-1 : 2012                                       | Labels, IFU                        | Section 4.2                    |
| <b>8.3</b>                            | Danger Symbol  | Yes       | EN ISO 15223-1 : 2012                                       | Labels                             | Section 4.2                    |
| <b>8.4</b>                            | Proper Label Format  | Yes       | EN ISO 15223-1 : 2012                                       | Labels                             | Section 4.2                    |
| <b>8.5</b>                            | Intended Purpose   | Yes       | EN 1041 : 2008  | IFU                                | Attachment# 2                  |
| <b>8.6</b>                            | Device and Component Identification                                  | Yes       | EN ISO 15223-1 : 2012                                       | Labels                             | Section 4.2                    |
| <b>8.7</b>                            | Instructions for use   | Yes       | EN 1041 : 2008<br>EN ISO 18113:2011                         | IFU                                | Attachment # 2                 |

**3. RISK ANALYSIS**

See Attachment 1: Non-Critical Disease Rapid Test Risk Management Report  
The Risk of this product is negligible and no further action needs to be taken.

**4. LABELS & INSTRUCTIONS FOR USE**

**4.1) Instructions for Use:** See Attachment 2: I.F.U.

**4.2) Direct Labeling of Product**

**See attached labeling**

### 5. MATERIALS SPECIFICATIONS

#### 5.1) Material Specification

| Materials                         |                            | Specification   |
|-----------------------------------|----------------------------|---|
| Test strip                        | Overall                    | 64 ± 0.5 mm x 3.5 ± 0.3 mm with five subcomponents  |
|                                   | a. Sample pad              | 31 ± 1.5 mm x 3.5 ± 0.3 mm<br>Absorb water within 1 second  |
|                                   | b. Conjugate pad           | 4-5 mm x 3.5 ± 0.3 mm<br>containing anti- <i>H. pylori</i> Antibody-colored particle conjugate  |
|                                   | c. Nitrocellulose membrane | 20 ± 0.5 mm x 3.5 ± 0.3 mm,<br>Flow rate: 4cm / 110 - 165 seconds<br>backing spotted with <i>H. pylori</i> Antibody (T line), and control reagent (C line).   |
|                                   | d. Absorbent pad           | 15 ± 0.7 mm x 3.5 ± 0.3 mm  |
|                                   | e. Vinyl matte adhesive    | 64 ± 3.3 mm x 3.5 ± 0.3 mm<br>Pass flow rate test at 45 °C for 7 days   |
| Anti- <i>H. pylori</i> Antibody 1 |                            | Purified antibody   |
| Anti- <i>H. pylori</i> Antibody 2 |                            | Purified antibody   |
| Colored particle                  |                            | Standard colored particles  |
| Conjugate pad fabric              |                            | Standard pad fabric   |
| Plastic Cassette                  |                            | 72 ± 4.1 mm x 20 ± 1 mm<br>A sample receiving well labeled S<br>The position of the T and C lines are marked on the cassette<br>Strip groove fits 3.5 mm x 64 mm strip<br>Pass flow rate test with specimens  |
| Pouch                             |                            | 120 ± 0.6 mm x 65 ± 3.2 mm size with 3-5 mm sealing margin<br>One side is printed with company information and test type.<br>The other side is blank for labeling of production information.<br>No dirty spots<br>Pass 3 day integrity of seal test |
| Desiccant                         |                            | 0.5 ± 0.1 g   |
| Package Insert                    |                            | Off white or white paper, 70g<br>A4 paper, left and right margin: 0.5 -1.0 cm; top and bottom margin: 1.0-1.5 cm.<br>Correct art work, color printing, no dirty spots   |

#### 5.2) Product Specifications

See Attachment 3: Product Specification

#### 5.3) Component Specifications

5.3.1) Test Characteristic

Nitrocellulose-based membrane strip with a T line pre-coated with unconjugated *H. pylori* antibody, a C line pre-coated with a control line antibody, and a conjugate pad containing colored particles conjugated anti-*H. pylori* antibody.

### 5.3.2) Kit composition and specifications

#### Kit Composition

| Composition                | Specification   |
|----------------------------|---|
| 1. Device                  | 25 single use devices in each kit. Each is individually sealed, and contains two items inside:<br>1. One cassette device composed of a test strip and a plastic housing cassette<br>2. One desiccant: 0.5 g |
| 2. Plastic dropper         | 10 plastic droppers in each kit for watery fecal specimens  |
| 3. Stool Collection Device | 25 devices in each kit<br>1 mL Extraction Buffer in each bottle   |
| 4. Package insert          | One insert each kit   |


#### Kit Box

| Specification                         | Dimension                               | Capacity        |
|---------------------------------------|---|-----------------|
| CTK H. Pylori Ag Specific (PM-R0192C) | 12.4 cm (W) x 22.2 cm (L) x 7.05 cm (H) | 25 test devices |

### 5.3.3) Test Appearance

The *Onsite* H. Pylori Ag Rapid Test is a cassette device. The device has following letters on the surface of the cassette:

- T: Test Line Position
- C: Control Line Position
- S: Sample Well



Both the Test line and the Control line in the result window are not visible before applying any samples.





If a positive specimen is applied, both the Test line and Control line will appear.



If a negative specimen is applied only the C line will appear. The Control line is used for procedural control. The Control line should always appear if the procedure is performed properly and the test reagents are working.



## 6. CLINICAL STUDIES & PERFORMANCE EVALUATIONS

See Attachment 4: Product Performance

## 7. STABILITY STUDIES

See Attachment 5: Stability Study

## 8. MANUFACTURING

### 8.1) Manufacturing Process

The entire manufacturing process has seven sequential process steps and is jointly accomplished by seven production groups

#### Step #1 Conjugate anti-*H. pylori* with colored particles and preparation of conjugation pad

- Preparation of colored particle solution.
- Add conjugation buffer, then add anti-*H. pylori* antibody to the solution
- Add blocker,
- Centrifuge
- Collect conjugate
- Dissolve the conjugate with conjugate suspension buffer
- Dispense the conjugates to the conjugate pad material
- Dry the conjugate pad


#### Step #2 Coat T and C line on the NC membrane

- Preparation of membrane lamination: Fix a 20 mm x 300 mm of nitrocellulose membrane onto the 300 mm x 66 mm of vinyl matte adhesive.
- Prepare the T line and C line coating solutions.
- Dispense the reagents to the T and C positions on the membrane with the coating machine (sprayer)
- The coated membrane is then dried

### Step #3 Lamination

- Assemble all the components to sheet according to scheme illustrated below
- Laminate the components according to scheme illustrated below
- QC samples the uncut sheet with the QC positive detection, specificity, and limit of detection panels.

[Cross Section Scheme of Lamination]



1. Vinyl matte adhesive
2. NC membrane coated with T and C line.
3. Conjugate pad
4. Absorbent pad

### Step #4 Cutting

- Laminated sheets are cut into the size of 3.5 mm x 68 mm.
- The size of the strip is inspected at the beginning, middle and end of the run.

### Step #5 Cassette Assembling and Sealing

- Pouch is labeled with name, mfg date, exp. date, and catalog, according to the documentation
- Assemble the cut strip into the plastic housing cassette and press to close the cassette
- Pack one cassette and one desiccant to each pouch
- Seal the pouch with the heat sealing machine
- Inspect packing process by reconciliation of the quantity of components picked, the quantity of pouch assembled, and the quantity of unused components.


### Step #6 Preparation of Sample Extraction Buffer

- Prepare Sample Extraction Buffer according to formulation
- Aliquot 1 mL into each extraction tube
- Ensure all bottles are properly labeled
- Inspect all the bottles for leakage at 2 Psi for 2 minutes.

### Step #7 Unit Packing

- Pick up all required components (Extraction Buffer, Devices) and labels including package inserts, zip lock bag (bulk package) and kit box
- Pack to unit
- Inspect packing process by reconciliation of the quantity of components picked, the quantity of units packed, and the quantity of unpacked components
- Packed product is moved to the quarantine area awaiting Final Acceptance Inspection

## 8.2) Manufacturing Flow Chart



### 8.3) Production Specifications

| <b>Production Process</b>               | <b>Specification</b>   |
|---|--|
| Pre-run Inspection                      | Antibody has less than 10% degradation by SDS-PAGE   |
| Colored conjugate solution              | Validated process  |
| Conjugation                             | First run inspection passes QC panel<br>Conjugate pad is uniformly sprayed with conjugates.<br>Dryness of the conjugate pad<br>Conjugate pad passes test with QC panel   |
| Coating                                 | The dispensing volume is calibrated at the start of the process by checking the accuracy of the volumes dispensed.<br>T and C lines are coated at designated area and inspected at the start run, middle run, and final run.<br>Dryness of the membrane<br>Coated membrane passes test with QC panel |
| Lamination                              | First run inspection passes test with QC panel<br>All the components are laminated correctly by visual inspection according to Inspection Sample Plan.   |
| Cutting                                 | 3.5 ± 0.1 mm wide each strip verified by inspection at start of run, middle-run, and final run.  |
| Cassette assembling and pouching        | All the strips are properly assembled into the cassettes by visual inspection during processing.<br>All the pouches are labeled correctly by visual inspection<br>All the components are packed correctly by reconciliation of the quantity of the components picked and used.                       |
| Preparation of Sample Extraction Buffer | Solution is made with the correct formulation.<br>Fill in volume 1 mL ± 0.1 mL based on the sampling plan<br>No leaking upon test inspection.  |
| Unit Packing                            | All the components are packed correctly by reconciliation of the quantity of the components picked and packed.   |

### 8.4) Documentation of Quality System

This product is manufactured in a facility certified to be in accordance with the ISO 13485 Quality System

See Attachment 6: ISO Certification

## 9.) **BATCH RELEASE CRITERIA**

Quality Control Inspection and Specifications:

Quality control inspection is performed on all raw materials and intermediate components produced, as well as the assembled kits. The inspection of the product is based on the following protocols:

### 9.1) **Incoming Material Inspection**

Incoming Material Inspections are performed on all in-coming raw materials. The quarantined raw materials are sampled per Inspection Sample Plan by QC department for inspection per specification described at section 5.1

Upon inspection, the QC supervisor will endorse the inspection data and assign quality status to the raw materials accordingly.

Only approved raw materials are used for production of product.

### 9.2) **In-Process Inspection**

All lots of individual intermediate components produced are inspected and tested during the process. A systematic sampling of each individual lot of intermediate components is taken for In-Process Inspection per specification described at section 5.1

Function tests are performed by assaying the intermediate components with the QC control panel, and with a reference intermediate component on the *OnSite H. Pylori Ag Rapid Test* if necessary. Physical inspections of the test components, such as the fill volume, appearance, and physical status are also carried out.

Upon analyzing the inspection results, the QC supervisor will assign the quality status of the intermediate components. The approved intermediate components are moved from the quarantined storage area to the approved storage area. Components that do not pass inspection will be rejected and not used in production.

### 9.3) **Final Acceptance Inspection**

Final acceptance inspection is carried out once all the components are assembled into the final packing unit. The inspection is to ensure only the product that meets the specification is released for distribution. The inspection includes:

- Document Inspection: Inspection of all production work records
- Physical Inspection: Inspect based on the Sampling Plan. The inspection includes checking the labels, lot number and expiration date of the individual components as well as the assembled kit. Inspection is also performed to ensure that all the components are packed.
- Performance Inspection: At least 40 tests submitted by the Production group are inspected with four *OnSite H. Pylori Ag Rapid Test* QC control panels. The panels consist of 6 members of a positive detection panel, 20 members of a specificity panel, 1 or 2 member of a precision panel and 3 members of limit of detection panel for testing product's positive detection, specificity, precision and limit of detection.

- Kits are ready for shipment once they have passed this final stage of QC inspection.
- QA will endorse the release of the product and retain at least 40 tests for future analysis.

#### 9.4) Procedure of the Final Acceptance Inspection

##### Introduction

Each kit of *Onsite H. Pylori Ag Rapid Test (R0192C)* contains the following components:

- a. 25 test devices, each sealed in a foil pouch with two items inside:
  - One cassette device
  - One desiccant
- b. 10 plastic droppers for watery fecal specimen
- c. 25 Stool Collection Devices
- d. One package insert (instruction for use)

Quality control evaluation is performed on every lot of these components.

##### Inspection Procedures:

The quality control procedures used for evaluation of the finished products are:

- a. Physical Inspection

According to the Inspection Sample Plan, obtain the required quantity of tests, inspect following parameters:

- Content in each package: Make sure each pack contains the correct quantity of components
- Labels: Make sure all labels correspond to documentation
- Pouch Integrity: Make sure the pouch is sealed properly
- Pouch Content: Open pouch to check contents. Make sure all the contents are included.
- Stool Collection Device Integrity: Make sure no leakage is visible

Record the number of the defects observed. Refer to the Acceptance Number (Ac). Pass inspection if the number (Ac) is less than the maximum allowable defects or defectives in a sample for the lot to be accepted based on the sampling plan.

- b. Performance Inspection

Each lot of the *Onsite H. Pylori Ag Rapid Test* is inspected for its positive detection, specificity, precision and limit of detection.

- Positive detection Inspection
  - Inspection is carried with the *OnSite H. Pylori Ag Rapid Test* Positive detection QC panel. The panel consists of 6 members, numbered SM-R0192-P1 to P11.
  - The assay is performed and interpreted using the procedure described in

the Product Insert.

- Each member specimen is assayed in duplicate.

- **Specificity Inspection**

- Inspection is carried with the *OnSite* H. Pylori Ag Rapid Test Specificity QC panel. The panel consists of 20 members, numbered SM-R0192-N1 to N20.

- The assay is performed and interpreted using the procedure described in the Product Insert

- **Precision Inspection**

- Inspection is carried out with *OnSite* H. Pylori Ag Rapid Test Precision panel consisting of one weak positive and one medium positive sample.

- The assay is performed and interpreted using the procedure described in the Product Insert.

- Each sample is assayed in 10 replicates. Total of 10 Devices are used for the inspection.

- The flow rate of each device is recorded during the inspection.

- **Limit of Detection (LOD)**

- Inspection is carried out with the 3 members of the LOD Panel.

- Each member specimen is assayed in duplicate.

- The assay is performed and interpreted using the procedure described in the Product Insert.

c. **Acceptance Criteria**

- The Positive detection Inspection result must meet the specification indicated in the following table.

| H. Pylori Ag Positive Detection Panel | Sample (N) | Result       |
|---------------------------------------|------------|--------------|
| SM-R0192-P1~P11                       | 11         | All Positive |

- The Specificity result must meet the specification indicated in the following table.

| H. Pylori Ag Specificity Panel | Sample (N) | Result        |
|--------------------------------|------------|---------------|
| SM-R0192-N1~N20                | 20         | 100% Negative |

- The Precision Inspection result must meet the specification indicated in the following table.

| H. Pylori Ag Precision Panel | Runs | Result                         | Flow Rate                |
|------------------------------|------|--------------------------------|--------------------------|
| SM-R0192-C- Weak Positive    | 10   | Equivalent test line intensity | migration: < 120 seconds |
| SM-R0192-C- Medium Positive  | 10   | Equivalent test                | migration:               |

|  |  |                |               |
|--|--|----------------|---------------|
|  |  | line intensity | < 120 seconds |
|--|--|----------------|---------------|

The Limit of Detection Inspection result must meet the specification indicated in the following table

| <b>H. Pylori Ag Limit of Detection Panel</b> | <b>Result</b> |
|--|---------------|
| SM-R0192-L1 (2 ng/mL)                        | +             |
| SM-R0192-L2 (1 ng/mL)                        | +             |
| SM-R0192-L3 (0.5 ng/mL)                      | + or –        |

Note: +: Positive; + or –, Positive or negative,

- d. Reference Components Used in QC  
The reference components used for inspection of individual intermediate components are approved components from previous production lots. These reference components are tested to ensure that results obtained are within the QC specification for the product.
- e. See Attachment 7: Certificate of Analysis.

**10.) CONCLUSION**

The *OnSite* H. Pylori Ag Rapid Test is developed, manufactured, and marketed according to the ISO13485 quality standard. In a study with specimens from 157 patient fecal samples tested on UBT gold standard method and the *OnSite* H. Pylori Ag Rapid Test had a Relative Sensitivity of 97%, a Relative Specificity of 94%, and an Overall Agreement of 95%. The test does not require equipment and can be performed by technicians with minimal training. It can be stored for 24 months at 2-30°C. Thus, the test is deemed acceptable for marketing and sale wherever the regulatory requirement is completed.