



FILTER PAPER FOR BREWERIES & BEVERAGE INDUSTRIES



Lasec® offers filters for malt analysis and malt filtration, according to **MEBAK, ASBC and EBC standards**, from world leading filtration and separation technology specialists Ahlstrom-Munksjö.

100/N FILTER PAPER
PRODUCT CODE: FLAS4328320



CHARACTERISTICS:

General purpose plain folded filter paper.

APPLICATIONS:

Filter paper for general filtration, malt analysis, malt filtration and degassing in beverage, brewery and sugar industries.

BENEFITS:

Low sodium, potassium and ammonium.

DIAMETER	BASIS WEIGHT G/M ²	SURFACE	FILTRATION SPEED S/10ML	HERZBERG S/100ML	TYPICAL RETENTION µM
320mm	85	Plain folded	30	120	6 to 8

E&OE

+27 21 531 7504 sales@lasec.com www.lasec.com

