


Lasec® introduces the latest “must-have” for your quality control lab:



THE WATERLAB WATER ACTIVITY METER



Take the guesswork out of food safety and shelf life testing with the NEW Waterlab Water Activity Meter by Steroglass from Lasec®.

This state-of-the-art Water Activity Meter directly measures the free water in products to determine the potential for microbial growth. Featuring a dew point sensor, large colour touchscreen, easy-to-navigate menu, and small footprint, the WaterLab is essential for any Quality Control lab.

Applicable industries include:

- ✓ Food Manufacturing
- ✓ Pharmaceutical and Cosmetics

Key features:

- Measuring range: 0.030 – 1.000 aW
- Accuracy: ± 0.003 aW
- Measurement time: < 5 minutes
- Certification: ISO compliant
21 CFR Part II compliant