



Accelerate your air monitoring process with the **AIRWEL® PLUS Air Sampler from Lasec®**

Many industries, including pharmaceuticals and healthcare, are subject to strict air quality regulations, such as GxP (Good Practice) and ISO standards which require routine environmental monitoring.

The AIRWEL® PLUS helps you with compliance by providing accurate and continuous air sampling data. Regular, automated microbial air sampling can prevent costly and time-consuming problems, such as contamination investigations, product recalls, or regulatory fines, in the future.



The AIRWEL® PLUS fitted with **(A)** a 90 mm plate and beside it an autoclavable stainless-steel grid and a wireless charging stand, **(B)** the disposable polyester grid head, **(C)** the optional tripod stand, and connected to **(D)** the gas station for compressed gas stream.

Contact your **Lasec® Sales Representative** today.
Ts & Cs apply.

+27 21 531 7504 sales@lasec.com www.lasec.com

SPECIAL FEATURES:

- ✓ **Compliant:** IP65, ISO14698, ISO17025, USP797, USP116, EU GMP/GLP.
- ✓ **Optimised Sampling:** Accelerated air sampling flow rates of 100 l / min or 200 l / min.
- ✓ **Adaptable:** Can be used with 90 mm or 55 mm Petri dishes.
- ✓ **Compact:** Weighing only 1.6 kg, it is light and easy to handle.
- ✓ **Portable:** The unit has a rechargeable battery which provides more than 10 hours of autonomy.
- ✓ **Hermetic and Hygienic:** IP 65 rated for dust and water and anti-germ painted.

EXCLUSIVE OPTIONS:

- ✓ **Auto-test:** Allows the user to check the airflow at any time.
- ✓ **Gas Station:** Allows the user to check the microbial quality of a gas source.



Contact your **Lasec® Sales Representative** today.
Ts & Cs apply.

+27 21 531 7504 sales@lasec.com www.lasec.com

 **Lasec®**

E&OE