

Aridity Science Laboratory Water Systems Overview

Let us help you find the perfect partner for your laboratory. Our portfolio encompasses systems that can provide Type 3 to Type 1 water, each with their own features. Use the helpful table below to narrow down the most suitable system for you based on your requirements.





MODEL	WATER QUALITY	PRODUCTION RATE (L/HR)	TOC MONITORING	RESERVOIR REQUIRED	FINAL FILTER OPTIONS	WALL MOUNTABLE	TOUCH SCREEN USER INTERFACE	REMOTE DISPENSER OPTION	DIMENSIONS (HXWXD)	DRY WEIGHT	KEY FEATURES	PAGE IN BROCHURE
Puro™ III	Type III	Up to 25, 45 or 85L/hr		✓		✓	✓		500 x 490 x 290mm 19.7" x 19.3" x 11.5"	23kg 50.7lbs	<ul style="list-style-type: none"> • RO dampening technology • Anti vibration • Variable flow rates 	10
Puro™ III T	Type III	10 or 20L/hr				✓	✓		500 x 485 x 330mm 19.7" x 19.1" x 13"	20kg 44.1lbs	<ul style="list-style-type: none"> • Compact • Integrated tank 	11
Puro™-i III	Type III	5, 10, 20 or 30L/hr	✓	✓			✓	✓	575 x 366 x 492mm 22.6" x 14.4" x 19.4"	23kg 50.7lbs	<ul style="list-style-type: none"> • Twin pass RO • Data capture • Drop by drop dispensing • Cartridge traceability 	12
Geno™ II	Type II	Up to 25 or 45L/hr		✓		✓	✓		500 x 490 x 290mm 19.7" x 19.3" x 11.5"	23kg 50.7lbs	<ul style="list-style-type: none"> • RO dampening technology • Anti vibration • Variable flow rates 	13
Geno™ II T	Type II	10 or 20L/hr				✓	✓		500 x 485 x 330mm 19.7" x 19.1" x 13"	20kg 44.1lbs	<ul style="list-style-type: none"> • Compact • Integrated tank 	14
Duo™ II I	Type II & I	Up to 25L/hr	✓	✓	✓	✓	✓	✓	500 x 490 x 290mm 19.7" x 19.3" x 11.5"	30kg 66.1lbs	<ul style="list-style-type: none"> • RO dampening technology • Variable flow rates • Drop by drop dispensing 	15
Duo™-i III I	Type III & I	5, 10, 20 or 30L/hr	✓	✓	✓		✓	✓	575 x 366 x 492mm 22.6" x 13.2" x 19.4"	23kg 50.7lbs	<ul style="list-style-type: none"> • Twin pass RO • Data capture • Drop by drop dispensing • Cartridge traceability 	16
Duo™-i II I	Type II & I	5, 10, 20 or 30L/hr	✓	✓	✓		✓	✓	575 x 366 x 492mm 22.6" x 14.4" x 19.4"	23kg 50.7lbs	<ul style="list-style-type: none"> • Twin pass RO • Data capture • Drop by drop dispensing • Cartridge traceability 	17
Alto™ I	Type I	Up to 2 L/min	✓		✓	✓	✓	✓	500 x 490 x 290mm 19.7" x 19.3" x 11.5"	23kg 50.7lbs	<ul style="list-style-type: none"> • Anti vibration • Automatic flush • Cartridge traceability • Drop by drop dispensing 	18
Alto™-i I	Type I	Up to 2 L/min	✓		✓		✓	✓	575 x 366 x 492mm 22.6" x 13.2" x 19.4"	23kg 50.7lbs	<ul style="list-style-type: none"> • Application specific cartridges • Data capture • Drop by drop dispensing • Cartridge traceability 	19

* Optional regional differences

* Dimensions for all accessories can be found on pages 20-21

SECTOR KEY

Look for the icon that matches your needs to find products designed with you in mind.

-  Academia/Government
-  Applied
-  Industrial
-  Pharma/Biotech