

WATER PURIFICATION



Lasec® is proud to be working with **Avidity Science** to meet your **Water Purification** needs.

Get your **ASTM type III, II** or **Ultrapure Water** using reverse osmosis, deionisation and ultrafiltration methods with the various Avidity Science water purifiers.



FEATURES AVAILABLE *

- Twin-pass RO to provide less water wastage
- RO sanitisation and optional tank UV lamp to prevent the risk of contamination and incorrect results
- Touchscreen display indicates water quality
- Dual wavelength ultrapure lamp

* *Dependent on model*

E-SERIES	I-SERIES
Easy to install, maintain and use	Innovative, high-level user experience
Economic and sustainable	Information on cartridge replacements with alarms and monitoring capabilities
Quality delivered	

Whatever your requirements, Lasec® and Avidity Science have a solution for you.

[Click here to complete our questionnaire](#)



Contact your **Lasec® Sales Representative** today.
Ts & Cs apply

+27 21 531 7504 sales@lasec.com www.lasec.com