



## BIOSAFETY CABINET **SERVICING & VALIDATION**

Lasec® boasts a team of service engineers led by an NSF accredited biosafety cabinet field certifier. NSF accreditation is viewed as the premier credential for individuals to certify Class II biosafety cabinets in hospitals, laboratories and research facilities.

Regular maintenance, validation and certification of a biological safety cabinet by a qualified service engineer guarantees that the cabinet will work as expected to protect you, your laboratory, and your surrounding environment from harmful exposure to biological agents.

### BIOSAFETY VALIDATION BY LASEC®

- Filter integrity leak test
- Airflow validation (inflow & downflow)
- Smoke test
- Visual inspection
- Levelling and integrity of base check
- Power installation check
- Functioning of display panel and controls test

**REGULAR SERVICING AND MAINTENANCE OF A BIOSAFETY CABINET IS ESSENTIAL TO PERSONNEL SAFETY**

**Lasec® - Trust our experience**

### **BIOSAFETY SERVICE ENGINEER TRAINING PROGRAMME**

First level biosafety training course for service engineers

- Includes training on airflow measurements, HEPA integrity, Leak tests, Smoke tests