

CARDIOVASCULAR SOLUTIONS

Lasec® has a comprehensive range of products for measuring and monitoring blood pressure and heart rate readings for the diagnosis of medical conditions.



PULSE OXIMETER

- Fingertip and handheld models available
- Large LED screen with up to four-direction display
- Blood oxygen saturation (SpO₂) and pulse rate monitoring
- Supports data transfer and data review

ELECTROCARDIOGRAM (ECG)

- High-speed automatic measurement, analysis, and diagnosis
- 12-lead ECG: simultaneously signals and prints the ECG waveforms
- High-resolution and precision thermal printing system
- Power supply: AC/DC, built-in rechargeable lithium battery
- Supplied with disposable electrodes (6x chest electrode and 4x limb electrodes)





BLOOD PRESSURE MONITOR

- Fully automatic desktop electronic sphygmomanometer
- Selectable modes for adult, paediatric and neonates
- High-definition colour LCD
- SpO₂ measuring function
- Supports data transfer and data review
- Option: fingertip probe

PATIENT MONITOR

- 8, 10 or 12 inch screen models
- Low power consumption
- Non-invasive and invasive parameter options
- Standard parameters:
 - ECG
 - RESP
 - TEMP
 - SpO₂
 - NIBP
 - PR
 - HR
- Optional parameters: ETCO₂ and 2IBP
- Optional features: touch screen, Wi-Fi and 5Ah lithium battery
- 3, 5 or 12-lead ECG, ST measurement and arrhythmia analysis



Contact your **Lasec**[®] Sales Representative for more information.

+27 21 531 7504 sales@lasec.com www.lasec.com

