



# Gel Electrophoresis Consumables



EDTA powder from Bio-Rad provides excellent chelation to assist with your electrophoresis needs.

**R 1600**

Product Expires November 2024

Product Code	Description	Unit of Measure
BBRD1610729	≥99% pure EDTA (disodium ethylenediaminetetraacetate dihydrate) powder	500 g

Prepare your buffers to a physiological pH range of 7.3 to 7.5 using Bio-Rad pure Tris base.



**R 2260**

Product Expires March 2025

Product Code	Description	Unit of Measure
BBRD1610719	≥99.8% pure Tris base (Tris[hydroxymethyl] aminomethane or Trimethamine)	1 kg

Contact your Lasec® Sales Representative today.

Prices valid while stocks last  
Prices exclude VAT & Delivery  
Ts & Cs apply

+27 21 531 7504    sales@lasec.com    www.lasec.com



## GelDoc Go Gel Imaging System

Product Code	Description
IBRD12009077	GelDoc Go Gel Imaging System with Image Lab Touch Software

The GelDoc Go Gel Imaging System gives you a benchtop imaging solution in a compact and evolved package. Acquire high-resolution, publication-quality images of both nucleic acid and protein gels. Advanced optical design for a large imaging area allows you to image up to four mini gels simultaneously and retain a compact size. Easily set up and controlled with built-in Image Lab Touch Software. Powerful enough to move your research forward, but small enough to fit almost anywhere.

## Trays for use with GelDoc Go



### UV/Stain-Free Sample Tray

One included with system

Product Code	Description
IBRD1708271	For imaging fluorescence blots, stain-free gels and blots, and gels stained with ethidium bromide, SYPRO Ruby, Oriole, or SYBR Stains, includes Gel Alignment Kit (IBRD12012160).

### White Sample Tray

Product Code	Description
IBRD12012165	For gels and blots stained with Coomassie Blue, copper, silver, or zinc stains, includes Gel Alignment Kit (IBRD12012160).

### Blue Sample Tray

Product Code	Description
IBRD12012160	For nucleic acid applications with GelGreen or any SYBR Stain, includes Gel Alignment Kit (IBRD12012190) and XcitaBlue goggles (IBRD1708185).

Contact your **Lasec® Sales Representative** today.

Prices valid while stocks last  
 Prices exclude VAT & Delivery  
 Ts & Cs apply

+27 21 531 7504    sales@lasec.com    www.lasec.com



E&OE