



New Universal Transport Medium (UTM-RT®) Tube Design for Roche Cobas® 5800/6800/8800 Platforms

Key Features:


- Same UTM-RT® medium
- No glass beads
- Compliant with package insert, no validations required
- Workflow improvement: Direct loading onto Cobas® racks
- Faster turnaround time
- Reduced risk of mislabelling
- Increased operator safety – no need to remove the swab from the tube



3ml Round Bottom UTM Tube without Glass Beads:

STOCK CODE	DESCRIPTION	UOM
HCPN3U090N	Without a Swab	PK 300
HCPN3U091N01	With Regular FLOQSwab®	PK 300
HCON3U92N01	With Progressive Minitip FLOQSwab®	PK 300

Workflow with Standard UTM



Workflow with new Copan UTM designed for Roche Cobas® platforms



Contact your **Lasec® Sales Representative** today.

Ts & Cs apply

+27 21 531 7504 sales@lasec.com www.lasec.com

