



A one-size-fits-all solution for **Specimen Collection** and **Transport**

eSwab® is a multipurpose solution designed for **anaerobic, aerobic, viral** and **Chlamydia sample collection and transport**. The unique Liquid Amies Elution medium allows for optimal sample recovery for both clinical and research settings.

The single device allows for a wide range of applications including traditional and automated bacterial culture, NAAT, Gram staining and rapid antigen assays.

Collection Sites:

- Wound
- Lesion
- Nasopharyngeal tracts
- Rectal
- Vaginal
- Eye
- Urethral

Fields of Application:

- Cutaneous infections
- Respiratory infections
- MRSA / MSSA
- CRE / VRE
- Sexually transmitted infections



Contact your **Lasec® Sales Representative** today.
Ts & Cs apply

+27 21 531 7504 sales@lasec.com www.lasec.com

