



EGG YOLK EMULSION ON SPECIAL!

Are you a microbiologist looking to detect *B. cereus*, *S. aureus*, *C. perfringens*, or *C. difficile* pathogens in your food samples?

EGG YOLK EMULSION 50%

Growth supplement for the culture of pathogens in microbial analysis of food, feed, water and environmental samples, intended for use with Baird-Parker Agar, *Bacillus cereus* MYP & Pemba, *Perfringens* Agar Base (TSC), and *Clostridium Difficile* Agar Base (Brazier's)

- Ready-to-use liquid emulsion - 100 mℓ vials
- 50% egg yolk / 50% distilled water
- Store in 2 – 8 °C; to be used within 24 hours of opening
- Non-hazardous - for in vitro laboratory use only

**SPECIAL
PRICE R230
VALID WHILE
STOCKS LAST**



Contact your **Lasec® Sales Representative** today.
Prices exclude VAT & Delivery. Ts & Cs apply.

+27 21 531 7504 sales@lasec.com www.lasec.com

 **Lasec®**

E&OE