

# Latex Powder-free Examination Gloves

Our Lasec® latex powder free examination gloves are manufactured using environmentally friendly energy in a low-carbon production facility.

They provide a shield in the laboratory and field.

Latex Material	Colour
Natural Rubber Latex	White



**Used for:** General purpose, hospital, laboratory, cleaning.

## KEY FEATURES:

- Size ranges: XS - XL
- High elasticity, flexibility, and high comfort
- 5-year shelf life
- Textured at fingertip
- Packaging: 100 / box
- Powder-free
- Complies to ASTM D3578 & EN455
- CE marked
- MDR 2017/745 Class I



## Nitrile Powder-free Examination Gloves



Our Lasec® nitrile powder free examination gloves are manufactured using environmentally friendly energy in a low-carbon production facility.

They offer the barrier protection you need.

Nitrile Material	Colour
Synthetic Nitrile <i>(best alternative for latex allergies)</i>	Blue



**Used for:** General purpose, hospital, laboratory, cleaning.

### KEY FEATURES:

- Size ranges: XS - XL
- Flexibility and high comfort
- 5-year shelf life
- Textured at fingertip
- For latex sensitive skin
- Powder-free
- Packaging: 100 / box
- Complies to ASTM D6319 & EN455
- CE marked
- MDR 2017/745 Class I



Contact your Lasec® Sales Representative for more information

+27 21 531 7504    sales@lasec.com    www.lasec.com



E&OE