

ZeptoMetrix®
an antylia scientific company

MOLECULAR CONTROLS AND PANELS FOR INFECTIOUS DISEASE TESTING

Infectious disease continues to be the **major cause of mortality in Africa**. Regular use of molecular controls aids in monitoring testing performance and enhances confidence in the quality and reliability of clinical laboratory operations.

Lasec® provides a variety of **ZeptoMetrix® NATtrol™ infectious disease molecular controls** for **routine quality control testing**. This extensive product offering provides easy-to-use whole organism controls designed to safely and efficiently assure that testing instruments provide accurate results.

NATtrol™ Infectious Disease Controls & Panels

- Epstein-Barr Virus (EBV)
- Herpes Simplex Virus (HSV)
- Human Hepatitis B & C Viruses (HBV & HCV)
- Human Immunodeficiency Virus (HIV)
- Meningitis / Encephalitis (ME)
- West Nile Virus (WNV)



ZeptoMetrix® NATtrol™

- ✓ **Convenient:** Ready-to-use
- ✓ **Cost-effective Shipping:** Utilising gel packs instead of dry ice for shipment
- ✓ **Innovative:** Whole organisms designed to mimic a clinical specimen
- ✓ **Efficient:** Full process controls to monitor extraction, amplification and detection
- ✓ **Safe:** Inactivated, and non-infectious
- ✓ **Versatile:** Can be used across various molecular testing platforms and assays

Contact your **Lasec® Sales Representative** for more information

+27 21 531 7504 sales@lasec.com www.lasec.com

 **Lasec®**

E&OE