

Every user of an autoclave has specific duties under the **Pressure Equipment Regulations (PER)** of the **Occupational Health and Safety (Act 85 of 1993)** to ensure the safety of the personnel using the equipment.

BASIC USER CHECKLIST FOR AUTOCLAVES



1

The autoclave must have a **Certificate of Manufacture**, countersigned by an **Approved Inspection Authority (AIA)** and **Importer's Declaration of Compliance**.



2

Equipment must be fitted with a **PER compliant nameplate** and **relevant safety accessories**.

*** The name plate must stipulate autoclave details including certification code and SANS 347 Hazard Category.**



3

Inspections with an **Approved Inspection Authority (AIA)** done **upon installation** with **re-testing** performed **every 36-months (3 years)** and reports kept on file.



4

Any **modification, repair, alteration records** must be received and **available for review** upon request by the **inspector**.



5

If your autoclave has **not been inspected and recertified** within the **last 3 years** – **contact Lasec® for immediate assistance**.

