



AFFORDABLE **GLASSWARE** TO MEET EVERYDAY LABORATORY NEEDS

Featuring **High quality, ISO Accredited, Borosilicate 3:3 glass (class A)**



CONTENTS

GENERAL LABORATORY GLASSWARE	02
VOLUMETRIC GLASSWARE	04
MICROSCOPY	06
FILTRATION GLASSWARE	07

GENERAL LABORATORY GLASSWARE

BEAKERS

- Short and tall form beakers with spout
- Uniform wall thickness
- Mechanically strong
- Easy to read white enamel printing and graduations
- Autoclavable at 121 °C



BOTTLES

Dropping

- Amber or clear
- Rubber teat or glass stopper

Media

- Amber or clear
- Round or square
- Blue polypropylene GL 45 screw cap and pouring ring

Reagent

- Amber or clear
- Narrow or wide neck
- Glass or stopper

Weighing

- With ground in stopper

CRUCIBLES

- Porcelain or Class A Borosilicate 3:3 glass
- Low or tall form
- Glazed options available
- Lids available



DISHES

Crystallising

- With or without spout

Petri

- Soda lime glass or premium borosilicate
- Lids included
- Available in volumes ranging from 20 - 3500 ml



FUNNELS

Filter

- Long or short stem

Buchner

- Glass - sintered disc with various porosities available

Separating

- Glass, teflon (PTFE) stopcock options
- Pear or conical shapes available
- With polypropylene stopper

MORTARS AND PESTLES

- Available in porcelain or glass
- Inside grinding surface of mortar and the rounded end of the pestle are unglazed for precise grinding and crushing.
- Mortar features a contoured spout for pouring and transferring of crushed material



TUBES

Culture

- Disposable
- With screw cap
- Ungraduated

Capillary

- Heparinised options available
- Sigillum wax available

Centrifuge

- Round or conical bottom
- Uniform wall thickness
- With or without screw caps
- With or without graduations

Test

- Soda lime or Class A Borosilicate 3:3 glass
- With or without rim
- Round bottom
- Graduated or ungraduated



WEIGHING SCOOPS

- Available sizes: 3 ml; 6 ml, 10 ml
- Made of borosilicate glass with open tubular arm, suitable for weighing small quantities of powder and dyes.



VOLUMETRIC GLASSWARE

BURETTES

Dr. Schilling

- Calibrated to contain (TC) the exact volume at the standard 20 °C reference temperature
- Automatic zero adjustment
- Complete with plastic bottle and stable plastic base

Straight Bore

- Calibrated to contain (TC) the exact volume at the standard 20 °C reference temperature
- Glass, teflon (PTFE) or rotaflo stopcock option
- Batch or individual certification available
- Graduations or printed in white or blue ink



CYLINDERS

Measuring

- Spouted or stoppered
- Hexagonal or round base
- Batch or individual certification

Nessler

- Blue, Ceramic Graduations - for Colorimetric Analysis
- Top rack dishwasher safe and autoclavable
- Available in 100 ml or 50 ml

FLASKS

Erlenmeyer

- Wide or narrow neck
- Plain, ground neck or screw cap

Filter

- With side arm
- Plain or ground neck

Iodine


- With ground neck and stopper

Kjeldahl

- Plain or ground neck

Volumetric

- Amber or clear
- Plastic or glass stopper
- Batch or individual certification available



PIPETTES

Pasteur

- Soda lime glass
- With or without cotton plug
- Sizes available: 150 mm or 230 mm

Serological

+

Volumetric

- Calibrated to contain (TC) the exact volume at the standard 20 °C reference temperature
- Batch or individual certification

The product images shown are for illustration purposes only and may not be an exact representation of the product.

E&OE

MICROSCOPY

COUNTING CHAMBERS



Fuchs-rosenthal

+

McMaster

+

Neubauer-improved

- With CE mark
- Approved for In-Vitro application according to IVD directive 98/79/EC

COVERSLIPS

Rectangular

+

Square

+

Round

- With CE mark
- Chemically resistant borosilicate glass
- Absolutely colourless, perfectly clear
- Suitable for fluorescence microscopy
- Thickness (0.13 mm - 0.16 mm)
- Approved for In-Vitro application according to IVD directive 98/79/EC
- Contains batch number for comprehensive information and traceability



SLIDES

Histobond

- Pre-cleaned
- Ground edges
- Autoclavable
- Permanent positive charge
- Dimensions: 76 x 26 mm
- Ideal for brain or spinal cord samples

Classic

- Pre-cleaned
- Boxes of 50
- Soda lime glass
- Compliance with DIN ISO8037/1
- Dimensions: 76 x 26 mm
- Thickness: 1 mm

Star frost

- Solvent resistant
- Colour coded marking area available in blue, green, pink, yellow and white



FILTRATION GLASSWARE



ALL GLASS VACUUM FILTER HOLDER WITH SILICONE STOPPER


- 300 ml borosilicate glass funnel
- Glass base & tubulated cap, 47 mm with silicone stopper
- Anodised aluminium spring clamp, 47 mm

* Please enquire about our range of filters

VACUUM FILTER HOLDER SET

- 300 ml Borosilicate glass funnel
- Glass base, 47 mm - fritted, fritted Teflon-coated or PTFE gasket
- Anodised aluminium spring clamp, 47 mm
- Silicone stopper

* Please enquire about our range of filters



**SCAN TO EXPLORE
OUR GLASSWARE
RANGE**



**CAPE TOWN
HEAD OFFICE**

52 Old Mill Road,
Ndabeni,
Cape Town,
7405

+27 21 531 7504

JOHANNESBURG

48 Gazelle Crescent,
Gazelle Place,
Corporate Park South,
Randjiespark,
Midrand,
1685

+27 11 052 2500

DURBAN

Unit 3, Windsor
Court, Ground Floor,
AFS House, 4 Derby
Place, Derby Downs,
Westville,
3629

+27 31 702 6001

PORT ELIZABETH

Unit 5.1, Building 3,
Ascot Office Park,
Block 7,
Conyngham Street,
Greenacres,
6001

+27 41 373 9173

BLOEMFONTEIN

Genius Loci
Building 5,
6 CP Hoogenhout
Street,
Langenhovenpark,
Bloemfontein, 9330

+27 51 447 9500



✉ sales@lasec.com

✉ international@lasec.com

🌐 www.lasec.com