

Lasec® HPLC Products


Is this happening in your laboratory?

- ▶ Eluents getting contaminated by dust and moisture in the laboratory air?
- ▶ Solvents evaporating resulting in fire and health hazards to staff?
- ▶ PFAS contaminating your tests?



HPLC Supply - Solvent Safety Caps

The Standard for Safe solvent supply in HPLC/UHPLC labs




Speak to your sales person about our useful Starter Kits suitable for most HPLC applications.

HPLC Disposal - Safety Waste Caps

Filter solvent vapours and ensure pressure equalization in the collection container.

✓ Resistant to Aggressive Chemicals

Through the use of pure PTFE and PE-HD Safety Waste Caps are resistant to organic solvents, acids and alkalis.



Avoiding PFAS without compromising on quality and safety: SCAT GreenLine Safety Caps

Innovative use of PEEK as a completely fluorine-free high performance plastic with excellent chemical and thermal properties to avoid PFAS in production of safety caps.

A broad range of safety caps and waste disposal solutions are available

Let us perform an HPLC safety check in your lab and assist you with solutions today.

