


# EFFICIENT DNA PURIFICATION USING NEXTFLEX SPRI TECHNOLOGY



Lasec® is excited to introduce the **NEXTFLEX® NGS Cleanup Beads** from Revvity – a powerful SPRI magnetic bead-based solution designed for efficient nucleic acid purification in NGS workflows.

- Efficiently **purify nucleic acids** for NGS library preparation
- Reversible magnetic bead binding for **optimal DNA recovery**
- Removes **unwanted components** (primers, salts, adapters)
- Precise **DNA size selection** capabilities
- Compatible with **manual and automated** workflows
- Drop-in replacement for most SPRI beads (**no protocol changes needed**)

## NEXTFLEX NGS Cleanup Beads



1. Capture DNA using magnetic beads
2. Isolate beads from impurities
3. Rinse magnetic beads with 70 % ethanol solution to eliminate impurities
4. Release DNA from magnetic beads
5. Transfer eluted DNA to a fresh plate

| PART NUMBER       | SIZE   |
|-------------------|--------|
| NOVA-BEADS-1 ML   | 1 ML   |
| NOVA-BEADS-5 ML   | 5 ML   |
| NOVA-BEADS-50 ML  | 50 ML  |
| NOVA-BEADS-500 ML | 500 ML |

Contact your **Lasec® Sales Representative** for more information

+27 21 531 7504    [sales@lasec.com](mailto:sales@lasec.com)    [www.lasec.com](http://www.lasec.com)



E&OE