



The World Health Organisation has recommended the use of **alcohol-based hand sanitiser with an alcohol content of at least 60%** as an **effective method of preventing the spread of COVID-19 through contact**. The high alcohol content in hand sanitisers means they are categorised as **flammable solvents**, so measures need to be taken to reduce the risk when storing bottles of hand sanitiser. **Lasec® offers a cabinet solution with 30 minutes of fire resistance, that can safely store up to 200 litres of hand sanitiser.**



SPECIFICATIONS:

- Offers **30 minutes of fire protection**
- Manufactured according to the **European Standard EN 14470-1** and the **South African standard SANS 54470-1** for fire safety storage cabinets

INTERIOR:

- 3 x shelves with a carrying load of 75 kgs
- 1 x bottom collecting sump and perforated insert

DIMENSIONS:

W x D x H (mm) 864 x 620 x 1947

COLOUR:

- Body colour - anthracite grey
- Wing door in warning yellow (RAL 1004)

E&OE

+27 21 531 7504 sales@lasec.com www.lasec.com