

Latex Powder-free Examination Gloves

Lasec® latex powder-free examination gloves are your trusted shield against **microbiological hazards** and **viral penetration**, ensuring reliable safety. Our gloves offer outstanding flexibility and stretch.

Latex Material	Colour
Natural Rubber Latex	White



Used for: General purpose, hospital, laboratory, cleaning.

The gloves are manufactured using environmentally friendly energy in a low-carbon production facility.

KEY FEATURES:

- Size ranges: XS - XL
- High elasticity, flexibility, and high comfort
- 5-year shelf life
- Textured at fingertip
- Packaging: 100 / box
- Powder-free
- Complies to ASTM D3578; EN455 & EN374-5:2016
- CE marked
- MDR 2017/745 Class I
- PPE Reg. 2016/425 CAT III



Contact your **Lasec® Sales Representative** for more information

+27 21 531 7504 sales@lasec.com www.lasec.com



Nitrile Powder-free Examination Gloves



Trust Lasec® nitrile gloves for enhanced protection against **chemicals**. Our nitrile gloves are the perfect alternative for those with latex allergies. Experience the flexible stretch and exceptional comfort of our powder-free examination gloves.

Nitrile Material	Colour
Synthetic Nitrile <i>(best alternative for latex allergies)</i>	Blue



Used for: General purpose, hospital, laboratory, cleaning.

The gloves are manufactured using environmentally friendly energy in a low-carbon production facility.

KEY FEATURES:

- Size ranges: XS - XL
- Flexibility and high comfort
- 5-year shelf life
- Textured at fingertip
- For latex sensitive skin
- Powder-free
- Packaging: 100 / box
- Complies to ASTM D6319; EN455 & EN374-5:2016
- CE marked
- MDR 2017/745 Class I
- PPE Reg. 2016/425 CAT III



Contact your **Lasec® Sales Representative** for more information

+27 21 531 7504 sales@lasec.com www.lasec.com

