



Introducing METAL-FREE TUBES

labcon®

In a world saturated with trace metals, certain applications require added assurance. Our Labcon® MetalFree® Centrifuge Tubes provide a pristine, metal-free processing environment which is crucial for experiments that are sensitive to trace metal contamination.

Crafted with precision for critical testing applications, Labcon® MetalFree® Tubes are not only free of ultra-low levels of trace metals, they are also lot tested and certified free of detectable Nuclease, Endotoxin, and ATP contamination.

These MetalFree® Tubes are a superior and reliable choice for the following samples:

- ✓ Environmental
- ✓ Food
- ✓ Chemical
- ✓ Biological


Contact your Lasec® Sales Representative for more information

+27 21 531 7504 sales@lasec.com www.lasec.com

 **Lasec**®

FEATURES:

- 100% medical-grade CFR21 compliant ultra-low metal virgin resins
- Certified below one part per billion (<1ppb) for trace metals
- Freezable down to -80°C
- Autoclavable
- Available in **1.5**, **15**, and **50 ml** sizes
- Choose from **resealable bags**, **cardboard racks**, or **bulk** packaging options



Elevate your lab processes to new heights with Lasec® and Labcon®.

Contact your Lasec® Sales Representative for more information

+27 21 531 7504 sales@lasec.com www.lasec.com



E&OE