

**Copan**  
innovating together

# Safely Collect, Transport, and Preserve Urine Samples



Copan® UriSponge™ is uniquely formulated to store urine samples for up to **48 hours** at room temperature, **reducing the risk** of overgrowth or overkill of specimens. It's unique design and formula introduces many advantages. The absence of boric acid increases operational efficiency by reducing the need for special handling during the transportation process whilst the sponge ensures that the correct volume of urine is absorbed, balancing the ratio between preservatives and urine.




## Key Benefits

- ✓ Boric Acid Free
- ✓ M40 – A2 Compliant
- ✓ IVDR
- ✓ 24 Month Shelf Life

Product Code	Description	Unit of Measure
HCPN802CE.A	UriSponge™ 12 x 80 mm	Pk 300

## How to use Practical and Reliable Collection



Copan Italia SPA. (2023). Product Focus: UriSponge™: Advanced active device for urine collection, transport, and preservation. Italy. Rev00